**METTLER TOLEDO**

## **InPro 2000 (i) pH Electrode**

### **Instruction Manual**

cs	Návod k obsluze	2
da	Brugsanvisning	6
de	Betriebsanleitung	10
en	Instruction Manual	14
es	Manual de instrucciones	18
fr	Instructions d'utilisation	22
fi	Käyttöohje	26
hu	Használati utasítás	30
it	Istruzioni per l'uso	34
ja	取扱説明書	38
ko	사용 매뉴얼	42
nl	Gebruikershandleiding	46
pl	Instrukcja obsługi	50
pt	Manual de instruções	54
ru	Инструкция по эксплуатации	58
sv	Bruksanvisning	62
th	คู่มือคำแนะนำ	66
zh	说明书	70

# InPro 2000 (i) pH elektroda

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### Obsah

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
## Instruction Manual

### Contents

1	Introduction .....	15
2	Safety instructions .....	15
3	Product description.....	15
4	Installation and commissioning .....	16
5	Operation.....	16
5.1	Calibration of the sensor and pH transmitter.....	16
5.2	Sterilization of sensor.....	16
5.3	Housing pressurizing.....	16
5.4	Temperature sensor .....	17
6	Maintenance .....	17
7	Trouble-shooting .....	17
8	Disposal.....	17
9	Guarantee.....	17

InPro and ISM are registered trademarks of METTLER TOLEDO Group in Switzerland, the USA, the European Union and a further five countries.

## 1 Introduction

METTLER TOLEDO pH sensors type InPro™ 2000 (i) are high-precision electrodes with liquid electrolyte, silver-ion trap and built-in temperature sensor. The InPro 2000 (i) sensors are suitable for the same wide range of applications as the corresponding 465 sensors without built-in temperature sensor, and provide an additional temperature signal for automatic temperature compensation of the pH signal during calibration and operation. InPro 2000 (i) sensors are in-situ sterilizable.

Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use.


## 2 Safety instructions

The sensor consists partly of glass and should be handled carefully. If the sensor is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn. Local safety regulations are applicable to the handling of the sensor.

We recommend that you only operate the electrode in combination with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

## 3 Product description

The inscription on each electrode contains the following information:

<b>METTLER TOLEDO</b>	Manufacturer of the electrode
<b>InPro™ 200Xi/SG/yyy/zzz</b>	Type designation; X=membrane glass; i=ISM Digital/SG=ancillary electrode/yyy=rod length in mm/zzz=type of temperature sensor
max. level	Maximum electrolyte level (to be respected when refilling the electrolyte)
Position electrode this side up	Position indication of sensor in the housing to avoid that electrolyte escapes through the refill opening
combination <b>pH</b>	Type of electrode
pH X...YY	pH measurement range
X...YYY°C	Temperature range for operation
Order No. 52 00X XXX	Order number
Refill XXXX	Electrolyte to be used for refilling (short number)
 SEV14 ATEX 0168 X	Identification and number of certificate
 IECEx SEV 14.0025X	Identification and number of certificate
 CML 22 UKEX 2108X	Identification and number of certificate
IS/I, II, III/1/ABCDEFG/T6	FM approved
See instruction manual!	Instruction (Observe the operating instructions)

In addition each electrode is allocated a serial number on the plug-in head to enable identification.

## 4 Installation and commissioning

1. On unpacking, check the sensor for possible mechanical damage of the pH sensitive glass, the diaphragm(s), the body and the connector. Report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the watering cap and rinse the sensor with deionized water. After rinsing, the sensor should only be dabbed dry with a tissue. Do not rub the pH sensitive glass, since this can lead to electrostatic charging and sluggish response times.
3. Carefully check the space behind the pH sensitive glass for the presence of any air bubbles and remove same by gently swinging the sensor in a vertical plane.
4. Remove completely from the electrode, the T-shaped rubber stopper in the electrolyte refill opening and its elastic securing band before inserting the electrode into the housing. Also remove possible electrolyte rests in the filling port.
5. Please refer to the instruction manual of the relative housing in order to install the sensor correctly.
6. Connect the sensor to the pH transmitter, using a corresponding cable. Please refer to the wiring and terminal diagram accompanying the cable.
7. For electrodes with Pt 100 temperature sensors, 3-wire technique is recommended in order to obtain best measurement accuracy. Please refer to the cable and transmitter instruction manuals for the corresponding cable wire connections.
8. Intelligent ISM® sensors such as InPro 200Xi enable "Plug and Measure" and extended diagnostics. For installation, commissioning and operation of the ISM system please observe the instructions for the transmitter, module and cable.

## 5 Operation

### 5.1 Calibration of the sensor and pH transmitter

For InPro 2000 (i) sensors, a 2-point calibration is recommended. The calibration can be made with the sensor either ready for mounting or already mounted within the housing. After removing the watering cap, the sensor is dipped in succession into two different buffer solutions with given pH reference values and the transmitter is calibrated to these buffer values.

### 5.2 Sterilization of sensor

In the case of sterile processes, such as in fermentation applications, the sensors are sterilized prior to the production cycle. Sterilization is normally effected by means of steam or superheated process medium with the sensor already mounted in the reactor inside its housing.

### 5.3 Housing pressurizing

During possible sensor sterilization and during normal operation, it is important that the housing and electrolyte pressure is adjusted to be higher than the effective pressure in the reactor vessel. This measure is necessary in order to avoid medium intrusion into the electrolyte. The optimal pressure difference value to keep the diaphragm clean and obtain acceptable electrolyte refill intervals has to be established empirically from case to case.

### 5.4 Temperature sensor

The integrated Pt 100/1000 temperature sensor is primarily intended for automatic temperature compensation of the pH signal during calibration and normal operation, and **not** for the accurate indication or control of the process temperature.

## 6 Maintenance

1. Regularly check the applied housing overpressure during operation, and re-adjust if required.
2. Regularly check the electrolyte level in the sensor. When reaching the end of the bulbous shaft section (reservoir), refill with the electrolyte indicated on the sensor up to the indication line for max. level.
3. In the case of electrodes with electrolyte bridge, reference electrolyte is first filled into the outer chamber (marked blue) up to the maximum level. The inner chamber (red marking) is then filled to a level approximately 0.5 to 1.0 cm higher than that of the outer chamber.
4. Carefully rinse the sensor tip and diaphragm(s) with deionized water after each production cycle. Under no circumstances must measuring solution be allowed to dry on these parts!
5. When the sensor is not in operation, store it with sensor tip and diaphragm well submerged into the same reference electrolyte as inside the sensor.
6. If the sensor is stored mounted in its housing, the same rule as in 4) above applies, but the storage electrolyte has to be slightly modified, with 2 parts of buffer solution pH 9.2 to 10 parts of electrolyte, this in order to avoid corrosion of housing parts.
7. If a sensor is stored dry for a few days by mistake, let it soak in the normal storage electrolyte for several hours before use.
8. Occasionally check the connector for possible traces of moisture. If necessary, clean well with deionized water or alcohol, and afterwards dry carefully.

## 7 Trouble-shooting

**At sluggish response time** (e.g. after a long and dry storage period), reactivation of the pH sensitive glass may be necessary. Dip the pH sensitive glass (but no further sensor parts) into Reactivation (HF) Solution (article No. 51 319 053) for 5 min. Afterwards, rinse well with deionized water and let sensor tip and diaphragm soak in reference electrolyte for minimum 12 hours.

**Rapidly increased resistance of reference electrode and/or drift** could be signs of protein precipitation. Use Cleaning Solution Pepsin/HCl (article No. 51 340 068) to remove such precipitation. At repeated problems of this type, consider switching to sensors with FriscoLyte™/9848 electrolyte, which counteracts protein precipitations.

## 8 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

## 9 Guarantee

On manufacturing defects, 12 months after delivery.

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  
 本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：  
 Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○ : 表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下  
 ○ : Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x : 表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。  
 x : Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

产品中有害物质的名称及含量 Toxic and hazardous substance name and containment in product						
部件名称 Part Name	有毒有害物质或元素 Toxic and hazardous substances					
	铅 Lead (Pb)	汞 Mercury (Hg)	镉 Cadmium (Cd)	六价铬 Hexavalent chromium (Cr6+)	多溴联苯 Polybrominated biphenyls (PBB)	多溴二苯醚 Polybrominated diphenyl ethers (PBDE)
插头 Elect. Connector	x	○	○	○	○	○
电路板 PCBA	x	○	○	○	○	○

本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：  
Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下  
○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。  
x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

For addresses of METTLER TOLEDO  
Market Organizations please go to:  
**[www.mt.com/contacts](http://www.mt.com/contacts)**

**METTLER TOLEDO Group**  
Process Analytics  
Local contact: [www.mt.com/contacts](http://www.mt.com/contacts)

Subject to technical changes  
© 08/2022 METTLER TOLEDO  
All rights reserved  
Printed in Switzerland. 52 001 519 C

Management System  
certified according to  
ISO 9001 / ISO 14001

**[www.mt.com/pro](http://www.mt.com/pro)**

\* 5 2 0 0 1 5 1 9 C \*