



Calibration and Service

OEM Parts & Specifications

Factory-Trained Technicians

Global Data Governance

Plans to Suit Any Need

Rainin 360+

Pipette Calibration

A Service Experience You Can Trust

METTLER TOLEDO

Pipette Performance Matters

Trust METTLER TOLEDO to Deliver

Your work is too important to risk inaccurate or inconsistent data due to poor pipette performance. Pipette service by METTLER TOLEDO is rigorous and comprehensive, executed by highly trained and experienced technicians who know your pipettes. Make certain these important research tools deliver the reliability and reproducibility your work depends on – with METTLER TOLEDO Service.

Reliable pipette service is part of the Rainin 360+™ system for optimal liquid handling results. Precision pipettes, high-purity tips and expert service are all backed by innovative tools, software and training to ensure exceptional performance across your workflows. See back cover to learn more.

With secure, global measurement standards and the world's largest service team, only METTLER TOLEDO offers the depth, breadth and level of process integrity and control your pipettes demand.



World's Largest Service Team

We have built the world's largest pipette service organization – including the largest network of ISO 17025-accredited service labs and onsite service teams – to ensure consistent, reliable service performance in over 40 countries.



Factory-Trained Technicians

Every pipette calibration technician is factory trained and equipped with state-of-the-art METTLER TOLEDO balances to assure accurate gravimetric measurements and provide a thorough and reliable pipette service experience.



Global Data Governance

Our proprietary Global Pipette Service Solution – GPSS – securely records measurement data directly from the balance to eliminate any potential for recording inaccurate results. Adhering to a single, global standard, it also links directly to PipetteX for electronic calibration certificate delivery. 24 CFR Part 11 compliant.

► Request pipette calibration today! www.mt.com/pipette-service



OEM Parts & Specifications

Maintenance and repairs are performed according to every major pipette manufacturer's specifications – from OEM parts to the recommended greases – to ensure that your instruments continue to perform with precision and accuracy.



Plans to Meet Your Needs

Whether you prefer onsite or mail-in service, we offer a wide range of service plans. With fast turnaround and the capacity to meet any laboratory's need – including ISO 8655 calibration service – METTLER TOLEDO has the capacity and expertise to meet your exacting needs.



Your Pipette Service Partner

We offer an experienced team of service professionals who can help you determine which calibration intervals and interim pipette verification tests best suit your needs and workflows.

Rainin 360+

Maximize the performance of your Rainin pipetting system with a platform to track pipettes and manage calibration, and verify the performance of any brand of pipette right at the bench. And make sure your whole team handles pipettes with proper technique – for any workflow – with pipette training from METTLER TOLEDO.



SmartCheck™

With Rainin SmartCheck, any user can quickly and easily verify that a pipette is performing correctly. With three test volumes using deionized water, SmartCheck delivers a pass/fail result in less than 60 seconds.

► www.mt.com/SmartCheck



PipetteX™ Asset Management Software

Download and store calibration certificates, and streamline day-to-day pipette tracking and calibration management. You'll spend a lot less time tracking and managing pipette calibration, and have more complete and accurate records. 21 CFR Part 11 compliant.

► www.mt.com/PipetteX



Good Pipetting Practice™

Improve your data quality with GPP™ – our comprehensive, systematic approach to maximizing pipetting accuracy and repeatability.

► www.mt.com/GPP

Rainin 360+, Global Pipette Service Solution, GPSS, PipetteX, SmartCheck, Good Pipetting Practice and GPP are trademarks of Mettler-Toledo Rainin, LLC.

www.mt.com/pipette-service

For more information

Mettler-Toledo Rainin, LLC
Laboratory Division
Local contact: www.mt.com/contacts

Subject to technical changes
© 07/2022 Mettler-Toledo Rainin, LLC. All rights reserved.

30805276 Rev A